

Unscrambler
ONLINE 

Installation Procedure
for
Unscrambler Online 2.2



Table of Contents

Installation Procedure for Unscrambler Online .2.2

1	Requirements prior to installation of Unscrambler-Online 2.2.....	1
2	How to install Unscrambler-Online 2.2	1
3	How to activate Unscrambler-Online 2.2	7

Table of Figures

Figure 1 - The setup.exe file in your installation folder	1
Figure 2 - Unscrambler-Online 2.2 Setup Wizard.....	2
Figure 3 - Unscrambler-Online 2.2 Setup Wizard.....	3
Figure 4 - Unscrambler-Online 2.2 Setup Wizard.....	4
Figure 5 - Unscrambler-Online 2.2 Setup Wizard.....	5
Figure 6 - Unscrambler-Online 2.2 Setup Wizard.....	6
Figure 7 - Shortcut to Unscrambler-Online in All Programs	7
Figure 8 - Unscrambler-Online 2.2 Activation Dialog.....	8
Figure 9 - Unscrambler-Online 2.2 Login Dialog	8
Figure 10 - Unscrambler -Online icon in Sys Tray.....	9

1 Requirements prior to installation of Unscrambler-Online 2.2

Prior to the installation of Unscrambler-Online, the user is requested to install .NET Framework 2.0. An installable version of .NET Framework 2.0 and its instruction manual are provided in a separate folder within the CD. Alternatively, .NET Framework 2.0 can be obtained from the following link.

<http://msdn2.microsoft.com/en-us/netframework/aa731542.aspx>

.NET Framework 2.0 can be installed on Windows 2000, 2003, XP, XP with SP1 & SP2 OS.

- Copy the contents of the folder Unscrambler-Online-Installer to any location, for example, <drive>:\Unscrambler-Online-Installer
- Move to the folder <drive>:\Unscrambler-Online-Installer
- Within the folder 'UO2.2' double click the 'setup.exe' file

2 How to install Unscrambler-Online 2.2

Step 1: Double click on the 'Setup.exe'




Figure 1 - The setup.exe file in your installation folder

Step 2: Double click on the ‘Setup.exe’




Figure 2 - Unscrambler-Online 2.2 Setup Wizard

Step 3: Select 'I Agree' and click on the 'Next'




Figure 3 - Unscrambler-Online 2.2 Setup Wizard

Step 4: Click on 'Next'. Default path for Installation is "C:\Program Files\Unscrambler-Online 2.2\"




Figure 4 - Unscrambler-Online 2.2 Setup Wizard

Step 5: Click on 'Next'




Figure 5 - Unscrambler-Online 2.2 Setup Wizard

Step 6: Click on 'Close'




Figure 6 - Unscrambler-Online 2.2 Setup Wizard

3 How to activate Unscrambler-Online 2.2

Step 1: Access the application through Programs >Unscrambler-Online 2.2 >Unscrambler-Online




Figure 7 - Shortcut to Unscrambler-Online in All Programs

Step 2: Provide the ‘Activation Key’ and Click ‘Activate’. If you do not have the Activation Key click on the “Obtain” button to receive a trial activation key to your email account that you have provided during product registration process. Paste the activation key and click on “Activate” button.




Figure 8 - Unscrambler-Online 2.2 Activation Dialog

Step 3: Provide ‘User Name’ and Click ‘OK’




Figure 9 - Unscrambler-Online 2.2 Login Dialog

