

UnscramblerX  
Predictor

# Installation Guide



UnscramblerX  
Predictor

# Unscrambler X Predictor Installation Guide

This guide is distributed with the media that includes an end user agreement, as well as the software described in it. Except as permitted by any such license, no part of this guide may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, recording, or otherwise, without the prior written permission of CAMO Software. Please note that the content in this guide is protected under copyright law even if it is not distributed with software that includes an end user license agreement.

The content of this guide is furnished for informational use only, is subject to change without notice, and should not be construed as a commitment by CAMO Software. CAMO Software assumes no responsibility or liability for any errors or inaccuracies that may appear in the informational content contained in this guide.

Version 10.1  
January 2011



©2011 CAMO Software. All rights reserved  
<http://support.camo.com/>

# Table of Contents

<b>1</b>	<b>Unscrambler X Predictor .....</b>	<b>2</b>
1.1	About the installation guide .....	2
1.2	System Requirements .....	2
1.3	Support.....	2
<b>2</b>	<b>Installing Unscrambler X Predictor.....</b>	<b>3</b>
2.1	Downloading the installer.....	3
2.2	End User License Agreement (EULA) .....	5
2.3	Enter serial number.....	6
2.4	Setting the Installation folder for Unscrambler X Predictor.....	7
2.4.1	<i>Installation progress</i> .....	8
2.4.2	<i>Completion of Installation</i> .....	9
<b>3</b>	<b>Activate the product.....</b>	<b>10</b>

## 1 Unscrambler X Predictor



The Unscrambler X Predictor is a CAMO initiative designed to make models generated off-line using The Unscrambler® X, more widely accessible to third party instruments and control systems. Unscrambler X Predictor can be interfaced into an existing software platform to provide users all of the predictive and diagnostic tools found in The Unscrambler® X, configured in a way that allows real time quality control.

### 1.1 About the installation guide

This manual, *Installation Guide*, describes the system requirements and explains how to install the software.

### 1.2 System Requirements

**[Caution: 64-bit Edition of OLUPX requires 64-bit processor]**

Windows 2003, XP (SP2 or higher), Vista, 2008, 7 or later

Processor: Pentium 233-megahertz (MHz) or higher

RAM: 1 GB minimum; higher is recommended

Hard disk space: 150 MB

.NET framework versions 2.0 to 3.5

Internet Explorer 7.0 or higher

### 1.3 Support

In case you experience difficulties in installing and activating the software license key visit <http://support.camo.com>

You may want to explore these other sources of support for The Unscrambler® X:

- [Learning Multivariate Data Analysis \(MVA\)](#)
- [Benefits to the use of MVA](#)
- [Webinars & Videos](#)
- [Web Workshops & Training Programs](#)
- [Suggested Reading](#)
- [Ask our Experts](#)

## 2 Installing Unscrambler X Predictor

### 2.1 Downloading the installer

Download Unscrambler X Predictor installer from <http://www.camo.com> and save it to a temporary folder.




Fig 1: Files and Folders in The Unscrambler® X Standalone folder

The windows installer wizard prompts the installation of Unscrambler X Predictor (Fig-2)





Fig-2 : The Windows Installer wizard

Once the installer loading is completed, a welcome dialog is shown where the user has the option to continue with the installation of the program or terminate the installation process (Fig-3)



*Fig-3 : Unscrambler X Predictor installation welcome dialog*

**Note:** Choose <Next> to proceed with the installation  
Choose <Cancel> to terminate the installation process.

## 2.2 End User License Agreement (EULA)




Fig-4: The End User License Agreement (EULA)


Users are requested to carefully review the EULA and convey the acceptance of the EULA to continue with the installation.

**Note:** Choose 'I accept' <Next> to proceed with the installation.  
If EULA is not acceptable, user can choose 'I do not' <Cancel> and abort the installation. Choose <Back> to step back to the screen. i.e. welcome dialog (Fig-3)

### IMPORTANT:

Choosing to agree indicates that you are “an authorized user and/or customer /end-user” and that you “have read the terms and conditions stated above” in the license agreement.

## 2.3 Enter serial number



Setup - OLUPX (64-bit)

**Personal Information**

Please specify your name and the company, then click Next.

Name:

Company:

Activation Key:

< Back   Next >   Cancel

Fig-5: The Enter serial number

## 2.4 Setting the Installation folder for Unscrambler X Predictor

Unscrambler X Predictor is installed to the location indicated in the text box. If you want to install the application to another location, type the new path or **<Browse>** to navigate to the new location.




Fig-6: Destination location for installation




Fig-7: Ready to install

## 2.4.1 Installation progress

Clicking on the **<Install>** button in the Ready to install dialog, the installation of the program begins with a progress dialog (Fig-8)




Fig-8: Installation progress indication

**Note:** Choose the **<Cancel>** button during the installation process to terminate the installation of the program.

## 2.4.2 Completion of Installation

Upon completion of installation process, a final dialog indicating the successful installation is displayed (Fig-9).




Fig-9: Completion of Installation

Choose the **<Finish>** button to complete the installation

### 3 Activate the product

Access the application through Programs >OLUPX (64-bit)10.1>OLUPX TestApp (Fig-9)




Fig-10: Unscrambler X Predictor location



[www.camo.com](http://www.camo.com)

### **CAMO Software**

NORWAY : Nedre Vollgate 8, N-0158, Oslo Tel: (+47) 223 963 00 Fax: (+47) 223 963 22

USA : One Woodbridge Center, Suite 319, Woodbridge, NJ 07095 Tel: (+1) 732 726 9200 Fax: (+1) 973 556 1229

INDIA : 14 & 15, Krishna Reddy Colony, Domlur Layout, Bangalore - 560 071 Tel: (+91) 80 4125 4242 Fax: (+91) 80 4125 4181

AUSTRALIA : PO Box 97, St Peters, NSW, 2044 Tel: (+61) 4 0888 2007

JAPAN : Shibuya 3-chome Square Bldg 2F, 3-5-16 Shibuya Shibuya-ku, Tokyo, 150-0002 Tel: (+81) 3 6868 7669 Fax: (+81) 3 6730 9539